

QR Code App

1. Download the FREE 3B Luna QR App from the App Store or Google Play to your smart device.
2. On the Luna II device, while in the **STAND-BY MODE**, briefly press the **ALARM SILENCE** button located to the left of the selector knob. On the G3 device, turn the knob until Report is highlighted in **BLUE**, push the knob, and select **iCODE QR/QR+**.
3. You will now see the iCode, iCode QR, and iCode QR+ screens. Select either iCode QR or iCode QR+ by turning the selector knob clockwise or counterclockwise. Once selected, the corresponding **QR Code** will appear in the display window. Using your smartphone or smart device, open the 3B Luna QR App. Once opened, select the iCode QR & iCode QR+ option.
4. Your device will now begin to **SCAN** for the QR Code. Center the QR Code within the square on your smart device.
5. Once scanned, you will see the Best 30 and **DETAIL** screens appear on your smart device. If applicable, the sleep data also was automatically transmitted and uploaded into the patient file associated with the device's Serial Number and Pin.
6. To exit from the iCode, iCode QR, and iCode QR+ screen on the Luna II, turn the selector knob to select **BACK** and press the knob; for G3 devices, press the **HOME** button twice.
7. For best results, use the QR+ option